



# STERILE PEDICLE SCREW SYSTEM:

Comprehensive offering of implants and specialized instrumentation for lumbar spinal fixation

- ✓ Reduction & Non-Reduction Screws
- ✓ Curved & Straight Rods
- ✓ Premium Instrumentation
- ✓ Compact Modular Trays



A variety of implant options are available, allowing flexibility based on the anatomy of each patient

**Screw Diameters Available:**

- 4.5 mm
- 5.5 mm
- 6.5 mm
- 7.5 mm

**Screw Lengths Available:**

- 20mm-60mm
- 5 mm increments



**Individually Packed**

All pedicle screws are individually sterile packed with one associated set screw



**100% Traceability**

Sterile packaging ensures product quality and traceability along with reduced hospital processing



**Zero Reprocessing**

Zero implant reprocessing helps eliminate the risk of infections

**About OC Spine**

OC Spine delivers the highest quality spinal implants and instruments that are derived from surgeon input. Headquartered in Savannah, GA, OC Spine develops innovative solutions produced with German manufacturing precision. OC Spine's goal is to provide sterile implant solutions with the objective of improving patient outcomes.

**Learn More**

 (888) 463-5803

 [Info@OCMedicalDevices.com](mailto:Info@OCMedicalDevices.com)